


# POLAR Walkthrough Private Bookmark

#### How to create a private bookmark?

1. Log into POLAR > Click Reports > Select any report you would like to access



2. Apply the required filters



3. Click on Bookmarks







### 4. Select Create Bookmark from Current Selections

mary PIP-QI Report Hel	p+	Patient List	Bookmarks <del>-</del>
	Pts with diagnosis DM without HbA1c recorded. Pts with intermediate risk of developing type 2 DM with high HbA1C Pts with intermediate risk of developing type 2 DM with high glucose who hi	as visited more than 10 times in the pa	ast 12 months.
nay be at risk of developing di t risk of developing diabetes (t	Create Bookmark from Current Selections		
cle of Care items. DRisk can be found via The Al	ustralian Type 2 Diabetes Risk Assessment Tool (AUSDRISK)		

## 5. Enter your bookmark title > Click on **Create**


- Clinical - 🖓	racked MRS Shared Health Summany DID-OI Report Help -	
	Create bookmark	X
At Risk Cohort	Title: Patient with active type 2 diabetes without HbA1c record in the last 12 months	
to help support t Risk Cohort" tab hknown) and thei	Description:	agement of t while the "M
e use of the AusE		isk Assessm
be available via	Close	to be indica
	HbA1c Manag Graph shows I	g <mark>ement</mark> atest HbA1c recor






#### How to access your private bookmark?

1. Log into **POLAR** > Click **Reports** > Select the report where the private bookmark was created



2. Click on **Bookmarks** > Select the bookmark you created from dropdown list



3. Related filters of the created bookmark will appear automatically



