

POLAR Walkthrough


QIPC – Patients with active chronic cardiovascular disease (CVD) with no total cholesterol results recorded in the last 12 months

Patient Cohort:

- Practice Active & RACGP Active (visited practice at least 3 times in 2 years)
- Active chronic CVD (Step 2)
- Total cholesterol not recorded in last 12 months (Step 3)


Note: Patient cohort can be modified by adding or removing steps according to your practice needs.

1. Log in to **POLAR** > click **Reports** > click **QIPC Clinic Report**



Practice Active and RACGP Active Filters are pre-applied.

2. Click **Clinical** > **CVD** > select **Active Chronic CVD**



3. Select the **Management** tab

Select date range **12 months** (NOTE: select date range before other filters)

On the Table, go to **Total Cholesterol** > select **No results** then ✓

The screenshot shows the POLAR QIPC Clinic Demo interface. The 'Management' tab is selected. The date range for pathology is set to '12 Months'. A table of patient data is displayed, with columns for various clinical metrics. The 'Total cholesterol' column is highlighted, and a dropdown menu is open, showing 'No results' selected. A red box highlights the 'No results' option in the dropdown.

Location	Chronic Disease Group	Most Seen Clinician	Usual GP	Gender	Age Group	Nicotine Replacement Therapy	Blood Pressure Lowering Therapy	BP	ACR	Lipid Lowering Therapy	Total cholesterol	HDL	LDL	Non-HDL	Triglycerides
Not Applicable	No	148/91	No results	Yes	3.3	1.3	No results	1.5							
Not Applicable	No	Not Recorded	No results	Yes	No	No results	No results	No results							
Not Applicable	No	Not Recorded	No results	Not Recorded	No	No results	No results	No results							
Not Applicable	No	Not Recorded	No results	Yes	No	No results	No results	No results							
Not Applicable	No	85/58	No results	Yes	2.7	1.2	No results	8.6							
Not Applicable	No	135/95	-0.6	Not Recorded	5.1	3.3	No results	8.39							
Not Applicable	No	173/88	No results	Not Recorded	3.9	1.7	No results	8.7							
Not Applicable	No	Not Recorded	No results	Not Recorded	2.2	1.2	No results	8.7							

4. Click on **Patient List** on the top right-hand side of the screen

The screenshot shows the top right-hand side of the screen. The 'Patient List' button is highlighted with a red box. Other buttons visible include 'Bookmarks', 'Selections', and a 'Patient Status' button.

5. Click on **Export to Excel** and save as an Excel file on your computer

The screenshot shows the POLAR interface. The 'Export to Excel' button is highlighted with a red box. Other buttons visible include 'Search', 'Filter', 'Print', and 'Share'.

For any further assistance in using POLAR, please contact your CQI officer or email cqisupport@swsphn.com.au